

Fast Refinery Gas Analysis on the Agilent 8890 Gas Chromatograph with Large Valve Oven

Author

Ian Eisele
Agilent Technologies, Inc.

Introduction

Refinery gases are mixtures of various hydrocarbons, permanent gases, and sulfur compounds. This complex mixture can be difficult to separate accurately and reliably. The Agilent 8890 GC fast refinery gas analyzer (RGA) provides a complete solution for the separation of refinery gases. Using a combination of packed, micropacked, and capillary columns, the analysis is separated into three channels for the simultaneous detection of hydrocarbons, permanent gases, and hydrogen with one injection (Figure 1). This configuration also benefits from the inclusion of the Agilent 8890 large valve oven. The large valve oven allows columns to be run isothermally, outside of the GC oven. This configuration improves oxygen response stability, as oxygen trends downward over time on porous polymers subjected to temperature programming.




Figure 1. 8890 fast RGA valve diagram. Columns 1, 2, and 3 are situated in the large valve oven.

The separation of hydrocarbons is performed on the FID channel with an Alumina PLOT column. Isomers of 22 hydrocarbons from C₁ to C₅ are separated, while C₆/C₆₊ hydrocarbons are backflushed from a DB-1 precolumn,

and elute as a group at the start of the chromatogram (Figure 2). Permanent gases are analyzed on the first TCD channel using a combination of HayeSep and MolSieve columns. These gases could include CO₂, CO, O₂, and N₂, with

the option of also analyzing H₂S and COS by adjusting the backflush time. The second TCD is dedicated to the analysis of hydrogen, running with nitrogen carrier gas for improved detection of hydrogen.




Figure 2. FID (A) and TCD (B,C) channels with compound identification.

www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2019
 Printed in the USA, March 1, 2019
 5994-0688EN