

DEC Caps for use with Gilson ASPEC™ Automated SPE Systems



- 1 and 3ml options
- Compatible with all Gilson ASPEC™ Systems
- Excellent sealing – needle and SPE Column
- Compatible with all 1 and 3ml Syringe Barrel SPE Columns
- Suitable for use with tab-less and conventional or winged/tabbed SPE columns

Introduction

Kinesis DEC Caps allow the automation of SPE methods using Gilson ASPEC™ Systems. Fully compatible with the majority of SPE columns, the Kinesis DEC Caps are inserted into the inlet of the syringe barrel style SPE column, allowing a positive pressure SPE method to be utilised.

1 and 3ml Options

Kinesis offers both 1 and 3ml DEC Caps. With the majority of analytical SPE methods using 1 or 3ml SPE Columns, the Kinesis DEC Cap supports most clinical, environmental, agrochemical and food solid phase extraction methods.

Tab-less & Conventional/Tabbed/Winged Columns

The outer diameter allows the Kinesis DEC Caps to be used with either conventional/tabbed/winged SPE Columns, as well as the tab-less version, often used in 96- and 48-well plate footprints.

Optimised Material

The proprietary material and design ensures an excellent seal with both the internal surface of the SPE Column and the needle of the appropriate Gilson ASPEC System.

Ordering Information

Part Number	Description	Pack Size
DC-GIL-001	Kinesis DEC Caps for Gilson ASPEC, 1ml	1000
DC-GIL-003	Kinesis DEC Caps for Gilson ASPEC, 3ml	1000

Samples

Samples of Kinesis 1 and 3ml DEC Caps are readily available. Please contact your local Kinesis office or distributor.



UNITED KINGDOM

Kinesis Ltd

Tel: +44 (0)1480 212122
 Fax: +44 (0)1480 212111
 sales@kinesis.co.uk
 Web: kinesis.co.uk

GERMANY, SWITZERLAND & AUSTRIA

Kinesis GmbH

Tel: +49 (0)2173 89 05-0
 Fax: +49 (0)2173 89 05-77
 sales@kinesisgmbh.de
 Web: kinesisgmbh.de

USA & CANADA

Kinesis Inc

Tel: (518) 289-5817
 Toll free in USA: (866) 934-6353
 Fax: (518) 289-5818
 sales@kinesis-usa.com
 Web: kinesis-usa.com

AUSTRALIA & NEW ZEALAND

Kinesis Australia Pty Ltd

Tel: +61 (0)7 3829 3996
 Fax: +61 (0)7 3829 3997
 sales@kinesis-australia.com.au
 Web: kinesis-australia.com.au

INDIA

Kinesis India Pvt. Ltd

Tel: +91-22-661 41 668
 Fax: +91-22-6668 1600
 info@kinesisindia.net
 Web: Kinesisindia.net

Kinesis recognise all trademarks and registered trademarks. E&OE



Kinesis Ltd kinesis.co.uk
Kinesis Inc kinesis-usa.com
Kinesis GmbH (formerly Abimed) kinesisgmbh.de
Kinesis Australia Pty Ltd kinesis-australia.com.au
Kinesis India Pvt. Ltd kinesisindia.net

 **Kinesis**
 Scientific Experts