



GERSTEL Solid Phase Extraction SPE

Vials and Caps

Description	Units	Part No.
10 mL Vial with screw neck	100	093640-038-00
10 mL amber vial with screw neck	100	093640-039-00
Screw cap, 18 mm, with septum silicone blue transparent, PTFE white, 35° Shore A, 1.3 mm	100	093640-040-00
1.5 mL Microliter vial, 32 x 11.6 mm, clear glass	100	093640-046-00
Magnetic screw cap, 9 mm, blue, with septum silicone white, PTFE red, 55° Shore A, 1.0 mm	100	015487-079-00
Cap for SPE waste position	1	015468-000-00
Septum for SPE waste position, 19 mm silicone blue transparent, PTFE white, 45° Shore A, 1.3 mm	10	093640-071-00

SPE Cartridges

Ready-to-use SPE cartridges fitted with transport adapters can be purchased from different suppliers:

- **Macherey-Nagel**
All M-N cartridges can be ordered directly
- **Varian**
Most standard cartridges can be ordered directly, other cartridges on request
- **Supelco**
Cartridges are available on request
- **LCtech**
Immunoaffinity cartridges are available on request



Syringes and syringe holders

Description	Units	Part No.
Disposable syringe needle for cartridge (Sterican G 27 x 3/4" 0.4 x 20)	100	009999-104-00
SPE syringe holder	1	015105-001-00
Microliter syringe 2.5 mL for SPE without needle	1	009980-063-00
GC needle, OD 0.63 mm, ID 0.36 mm, needle length 60 mm	2	009980-066-00
LC needle, OD 0.72 mm, ID 0.37 mm, needle length 60 mm	2	009980-067-00

SPE Station replacement parts and tools

Description	Units	Part No.
SPE Starter Kit	1	015547-000-00
Insert for 1 mL SPE cartridge for SPE Station (Adapter)	1	015279-001-00
Insert for 3 mL SPE cartridge for SPE Station (Adapter)	1	015279-003-00
Dismantling tool for cartridge inserts	1	015585-000-00
Insert for 2 mL Vials, SPE Station (Adapter)	1	015487-002-00
Insert for 4 mL Vials, SPE Station (Adapter)	1	015487-004-00
Insert for waste position for SPE Station	1	015212-000-00
Cap for SPE waste position	1	015468-000-00
Septum for SPE waste position, 19 mm silicone blue transparent, PTFE white, 45° Shore A, 1.3 mm	10	093640-071-00
Bleeder needle for SPE Station	1	015238-001-00
Silicone tube ID 6 x 1.5 mm for waste position	1	015994-002-00
Silicone tube ID 2 x 1 mm for needle vent	1	015995-002-00
Threaded bottle 1.5 L HDPE wide-necked, 108 x 108 mm, height 209 mm, thread 80 mm	1	015504-015-00
Screw cap PP for bottle with neck size 80 mm	1	015504-016-00
Service kit for SPE carriage	1	016074-001-00
Replacement glass balls for SPE carriage	10	016067-010-00
Replacement roller kit for SPE carriage	3	016059-000-00



GERSTEL certified SUPPLIES

SPE cartridge tray consumables

Description	Units	Part No.
Septa strips for SPE cartridge magazine with 72 positions (6 units used for one tray)	6	015590-006-00
Septa strips for SPE cartridge magazine with 98 positions (7 units used for one tray)	7	015590-007-00
Septum strip base for SPE cartridge magazine with 72 positions	1	015276-106-00
Septum strip base for SPE cartridge magazine with 98 positions	1	015276-101-00



GLOBAL ANALYTICAL SOLUTIONS

GERSTEL GmbH & Co. KG
Eberhard-Gerstel-Platz 1
45473 Mülheim an der Ruhr
Germany

GERSTEL, Inc.
701 Digital Drive
Suite J
Linthicum, MD 21090
USA

GERSTEL AG
Wassergrabe 27
6210 Sursee
Switzerland

GERSTEL K.K.
2-13-18 Nakane, Meguro-ku
152-0031 Tokyo
Dai-Hyaku Seimei Toritsudai
Ekimae Bldg 2F
Japan



+49 208 - 7 65 03-0
+49 208 - 7 65 03 33
@ gerstel@gerstel.com
www.gerstel.com

+1 410 - 247 5885
+1 410 - 247 5887
@ sales@gerstelus.com
www.gerstelus.com

+41 41 - 9 21 97 23
+41 41 - 9 21 97 25
@ gerstel@ch.gerstel.com
www.gerstel.de

+81 3 57 31 53 21
+81 3 57 31 53 22
@ info@gerstel.co.jp
www.gerstel.co.jp



Subject to change. GERSTEL®, GRAPHPACK® and TWISTER® are registered trademarks of GERSTEL GmbH & Co. KG.
© Copyright by GERSTEL GmbH & Co. KG