

Agilent's Industry-Leading GC/MS Systems

The right choice for reproducible results







Tools, Techniques, and Turnaround Times

You can't go wrong with a GC/MS instrument from Agilent. We start with the world's most trusted gas chromatographs. Our latest models are super smart: they tell you when to change a column, let you know if there's a leak, and collect system information to help you be more productive. Then we pair those GCs with quadrupole, triple quadrupole, and time-of-flight mass spectrometers—so you can choose the technique that best suits your needs. Whether you're doing routine testing or cutting-edge research, we have a system for that. And we have lots of ways to accelerate your turnaround time.



Agilent Intuvo 9000 GC

Systems



Agilent 5977B GC/MSD system

Agilent 5977B GC/MSD system

The source of new possibilities

Put the industry's most popular and powerful single quadrupole GC/MSD system to work for your lab—and watch productivity rise. This system is ideal for routine applications such as environmental, chemical, petrochemical, food, forensic, pharmaceutical, and material testing.



Agilent 7000D triple quadrupole GC/MS system

Agilent 7000D triple quadrupole GC/MS system

The definitive answer

The 7000D triple quadrupole GC/MS is the latest model of the most successful GC/MS/MS in history. The right choice for most applications, the 7000D is the fourth-generation version of the GC/TQ used by agencies around the world to create and validate many now-standard GC/MS/MS methods.



Agilent 7250 GC/Q-TOF system

Agilent 7250 GC/Q-TOF system

See the whole picture

Identify unknown compounds and uncover additional data using the unique capabilities of our quadrupole time-of-flight GC/MS system, with its low energy capable electron ionization (EI) source. This revolutionary all-in-one system is ideal for performing complex metabolomics studies, screening for pesticides in challenging matrices, identifying analytes in herbal extracts, monitoring contaminant levels in fuel, and more.




Agilent 7010B triple quadrupole GC/MS system

Agilent 7010B triple quadrupole GC/MS system

Expand your capabilities with sensitive, robust MS/MS

Meet the most sensitive version of our compact, benchtop triple quadrupole GC/MS. The 7010B triple quadrupole GC/MS delivers attogram detection limits in Elw mode. What's more, its breakthrough sensitivity lets you optimize sample preparation, reduce maintenance cycles by injecting less, and lower your detection limits.




Make Your Mass Spectrometry More Powerful

The Agilent MassHunter suite

Simplify data analysis, review, and reporting

Looking for an easy way to control the latest GC technologies—including sampling techniques and dedicated GC detectors? MassHunter Acquisition enables any user to produce great results with Agilent GC/MS systems.

Time-saving tools such as autotune, MRM Optimizer, and spreadsheet sequence import unlock the power of your instrument for method development and high-throughput workflows. What's more, optional compliance features provide peace of mind for labs that require the highest levels of data integrity.





Innovative Options

Agilent JetClean self-cleaning ion source

Clean less, analyze more

During routine analysis, matrix deposits inevitably build up. As a result, users have had to routinely remove the ion source, scrub the lens and other components, then put it all back together and recalibrate the instrument. We offer a better option. Using a carefully controlled hydrogen flow, the Agilent JetClean self-cleaning ion source automatically removes even the most stubborn deposits.




Agilent QuickProbe technology

Analyze samples in less than 60 seconds

With Agilent QuickProbe technology, your lab can quickly and directly analyze the composition of tablets, powders, and liquids to determine whether controlled substances are present. In fact, this real-time GC/MS technique enables you to identify compounds with little or no sample preparation.



Agilent CrossLab services



Learn more about Agilent CrossLab, and see examples of insight that leads to great outcomes, at www.agilent.com/crosslab

Learn more:

www.agilent.com/chem/gc-ms

Find a local Agilent customer center:

www.agilent.com/chem/contactus

U.S. and Canada 1-800-227-9770 agilent_inquiries@agilent.com

Europe

info_agilent@agilent.com

Asia Pacific

inquiry_lsca@agilent.com



