



Detect Antigen-Specific CD8⁺ and CD4⁺ T Cells with Confidence

Immudex MHC Dextramer[®] reagents are designed to obtain high-quality results. They contain multiple binding sites and fluorophores to boost avidity in binding antigen-specific T cells even with rare or low-affinity receptors.



MHC I Dextramer[®] reagents can identify low-affinity antigen-specific CD8⁺ T cells not detectable with the corresponding tetramer. Adapted from Dolton et al., 2014.

Try Our Special Offer!

MHC I and MHC II Dextramer[®] reagents are now available with a **10% discount**

To order your MHC Dextramer[®] please send an e-mail to

bio@perlan.com.pl by specifying your

- Cat. Number - Scan the QR code or visit:
www.immudex.com/mhc1 (for MHC I)
www.immudex.com/mhc2 (for MHC II)
- Fluorophore (PE, APC or FITC) and test size (50 or 150)
- Include your discount code: **IMM.2022**

Terms and Conditions

- Promotion is valid only from May to June in Poland
- Promotion is valid only for the products available on our website and test size reported on this flyer. Custom products are not included in this promotion.
- Contact us: **bio@perlan.com.pl**

Local distributor:



Perlan Technologies Polska Sp. z o.o.
ul. Puławska 303
02-785 Warszawa | Polska
bio@perlan.com.pl




MHC I Alleles



MHC II Alleles

Immudex U-Load Dextramer® - Load Your Own Peptide



U-Load Dextramer® reagents allow you to create custom Dextramer® from the comfort of your lab by loading **your own peptide** of interest on MHC I* and MHC II** receptive monomers.

* Immudex is the distributor of immunAware's easYmer® peptide-receptive MHC I monomers.

** Immudex peptide-receptive U-Load MHC II monomers.

Try Our Special Offer!

MHC I and MHC II U-Load Dextramer® Kit are now available with a **10% discount**

To order your U-Load Dextramer® Kit, please send an e-mail to

bio@perlan.com.pl by specifying your choice of MHC allele

- discover our extensive MHC allele list by scanning the QR code or visiting our website: www.immudex.com/u-load-alleles
- Test-Size (50, 150 or 500 for MHC I; 50, 150 or 500 for MHC II)
- Fluorophore (PE, APC, or FITC)
- Include your discount code: **IMM.2022**



Terms and Conditions

- Promotion valid only in Poland from May to June
- Promotion is valid only for the products available on our website and test size reported on this flyer. Custom products are not included in this promotion.
- Contact us: **bio@perlan.com.pl**

Local distributor:



Perlan Technologies Polska Sp. z o.o.
ul. Puławska 303
02-785 Warszawa | Polska
bio@perlan.com.pl

© Immudex ApS Denmark, 2022
For research use only. Not for use in diagnostic or therapeutic procedures.